HS Troponin Guidelines

		HS Trop Once	
Providers can order any single order or combination	HS Trop Series: B/L		
	HS Trop Series: 2hr		
	HS Trop Series: 6hr		
Collection Requirements	2hr	1 - 3 hours post B/L or Once collection	
Collection Requirements	6hr	5 - 7 hours post B/L or Once collection	

Issues lab staff can address without clinician input	A)	Re-timing HS Trop 2hr
	B)	Re-timing HS Trop 6hr

Timing Guidelines					
Timing problem	What to do	Documentation	Example		
B/L <i>more than 60 min</i> early or late	 Re-time 2hr to 2 hours post B/L collection AND Re-time 6hr to 6 hours post B/L collection 	Edit spec time (and date if necessary) in Meditech AND Replace original labels on board with new correct labels	Original orders B/L: 0000 2HR: 0200 6HR: 0600 B/L collected late at 0200 Retime: 2HR: 0400 6HR: 0800		
2hr more than 3 but less than 5 hours after B/L	Receive and perform 2hr *DO NOT RETIME 6HR *	Bench tech will add comment to delta result			
2hr 5 or more hours after B/L	 Cancel 2hr with comment AND Receive specimen as 6hr 	Cancel comment: "no specimen received"Explain delay if possible	"No 2hr specimen received due to pt in cath lab until 1700."		
6hr more than 7 hours after B/L	Receive and perform 6hr	Bench tech will add comment to delta result			

Clinical Decisions CANNOT be made by lab staff				
MUST record RN/provider Meditech ID or name in order/cancel comment	A)	Retiming HS Trop Once or B/L		
	B)	Cancelling HS Trop orders, including duplicates		
	C	Backordering for Inpatients		
	D)	Backordering for ED when blood in lab not recently collected		
	E)	Changing TYPE of Trop Order (ex: changing B/L to Once)		
See HS Troponin Order Issues Flow Chart for help				

HS Troponin Order Issues

